



Endo-L (Instruction For Use)

Please read carefully before use

INTENDED USE:

- For use as root canal lubricant during mechanical instrumentation in endodontic treatment of teeth.
- For use as a chelating agent during cleaning and shaping of canals in endodontic treatment of teeth.
- For use as a root canal conditioner during chemico-mechanical preparation of calcified root canal and shaping the canal, aiding in the softening of the calcifications and removal of debris from the prepared tooth.

COMPOSITION:

Ethylenediaminetetraacetic Acid (EDTA), Water.

PRODUCT DESCRIPTION:

- Endo-L is a ready to use EDTA Liquid, packed in Multi-dose plastic dispensing bottle. It is both lubricating and chelating material in liquid form, to be used during cleaning and shaping of the root canals during endodontic treatment. This is non sterile device & does not require sterilization while being used.
- Endo-L is a water based EDTA liquid specifically formulated to provide cleansing of the root canal preparation. Endo-L facilitates the shaping of the root canal while rotary and hand instrumentation.
- Endo-L is a plain EDTA liquid formulated to chelate the calcified surfaces in root canal and aid in shaping the canals by softening the calcifications.
- Endo-L is designed to be used with endodontic irrigation solutions such as sodium hypochlorite. This action allows for pulp tissue, dental shavings and debris to be easily removed from the root canal.
- Endo-L aids in achieving an appropriate internal canal surface for penetration of other medicaments used in endodontic therapy by cleaning and conditioning the root canal.

PROCEDURE & PRODUCT HANDLING:

- Open the pulp chamber in the usual manner and obtain adequate access to the root canal or canals.
- Extricate the pulp and flush canal with sodium hypochlorite solution.
- Dispense Endo-L either directly from the multi dose plastic dispensing bottle, using a syringe with irrigation tip or needle in and around the canal orifice, OR dispense Endo-L from the multi dose bottle, into a dappen dish and carry it to the canal using an irrigating syringe.
- The lubricant action is achieved on mechanical instrumentation, allowing the files and reamers to act smoothly. EDTA acts by chelating the calcium in the dentinal tubules there by aiding in the removal of the hard tissue.
- Using a small hand file, explore the canal, gently irrigate and evacuate with irrigating solution, thoroughly flush the canal with normal saline and irrigate with Endo-L.
- Apply Endo-L to the hand / rotary file or reamer and use it to desired depth. Repeat the instrumentation to achieve the desired canal shape and complete debridement.
- Once shaping is completed, thoroughly flush and evacuate the canal using Endo-L to remove all irrigants and smear layer.
- Dry the canal and obturate.

CAUTION / WARNING:

- Device to be used by a qualified dentist or dental technician only.
- Do not inject material under pressure into the canal.

PRECAUTIONS:

- Avoid contact of the device with skin and mucosa, avoid contact with eyes. Avoid inhalation of fumes. Do not take internally.
- Eye contact: If the product comes in contact with eye, it may cause irritation. In case of contact with eyes, rinse immediately with plenty of water and seek medical attention.
- Skin Contact: Avoid contact with skin, it may cause skin sensitization. If contact with skin occurs immediately wipe off thoroughly with cotton and then wash well with soap and water. If skin rash and sensitization or other allergic reaction occurs discontinue use and seek medical attention.

CONTRAINDICATIONS:

- Do not use the product if the patient is allergic to any of the ingredients.
- Do not use the product on mucosa or skin.
- Not Intended for use beyond the confinement of root canals of teeth.

INTERACTIONS WITH OTHER DENTAL MATERIALS:

No such Reactions Noted. Kindly communicate and inform manufacturer in details of any instance.

ADVERSE REACTIONS:

No such Reactions Noted. Kindly communicate and inform manufacturer in details of any instance.

PACKAGING:

Multi dose plastic dispensing bottle.

STORAGE CONDITIONS:

Keep Away from Direct Sunlight. Keep locked up and out of the reach of children.

Store at a temperature between 10°C to 30°C.

Shelf life (2 years from the date of manufacturing) not affected after opening, if stored in a cool dry place.

DISPOSAL CONDITIONS:

Dispose in accordance with local, state and federal environmental regulations.

This product is specially formulated for use in dentistry and to be used by qualified dental personnel.

Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

B-18, Dwarka Indl. Est., Naik Pada, Waliv, Vasai (East),
Dist. Palghar - 401 208, Maharashtra, India

www.maarcidental.com | care@shivaproducts.net | T: +919967806543

EC REP

Obelis S.A
Bd GénéralWahis, 53,
B-1030 Brussels, Belgium
Tel: +3227325954 | Fax: +3227326003
Email: mail@obelis.net

CE
0123

Revision	05
Date	30-11-2020

These instructions for use should be kept for the duration of product use.



Symbol for "BATCH CODE"



Symbol for "USE BY"



Symbol for "NON STERILE"



Symbol for "KEEP DRY"



Symbol for "CAUTION"



Symbol for "CONSULT INSTRUCTION FOR USE"



Symbol for "CATALOGUE NUMBER"



Symbol for "EUROPEAN COMMUNITY AUTHORIZED REPRESENTATIVE"



Symbol for "DO NOT USE IF PACKAGE IS DAMAGE"



Symbol for "KEEP OUT OF THE REACH OF CHILDREN"



Symbol for "MANUFACTURER NAME & ADDRESS"



Symbol for "TEMPERATURE LIMITATION"



Symbol for "KEEP AWAY FROM SUNLIGHT"

The material has been developed solely for professional dental use. Application should be carried out strictly according to the instruction for use. Liability cannot be accepted for damages resulting from failure to observe the instructions of the stipulated area of applications. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in this instruction sheet. Description and data constitute no warranty of attributes and are not binding.