


Blitz White



Known as: In-office Tooth Whitening kit

Description: It is a hydrogen peroxide based In-Office Tooth Whitening kit requiring minimal chair time. It contains desensitizers to maximize patient comfort.

Packaging:

- 2 x 0.4 g Powder
- 2 x 3 ml Liquid
- 1 x 3 g Gingiva Shield
- 2 x Micro Brushes
- 2 x Dispensing Tips

Product Code: 8801 / 001

Properties:

- In-office system with De-sensitizing agent
- Less post procedural sensitivity
- Quicker Whitening Result
- Contains Fluoride for Re-mineralization

Shelf life and storage:

- Best before 24 months under proper storage
- Store at a temperature between 10°C to 30°C.



Manufactured by:
Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com