

Calx-O



Known as: Calcium Hydroxide Paste (With Barium Sulphate)

Description: It is premixed high pH Calcium Hydroxide Water based paste with stabilizers and Barium Sulphate. To be used as a root canal dressing and filling material during endodontic treatment of teeth and as a temporary filling material for pulp capping and pulpotomies.

Packing:

1.5 g x 2 Syr with 3 Applicator tips
3 g x 1 Syr + 3 Applicator Tips (For Sale in India Only)

Product Code:

9102/1.5x2
9105/003 (For Sale in India Only)

Indications:

- Flowing Canal
- Apexification - Hard tissue barriers
- Root Resorption

Advantages:

Radio-opaque with excellent biocompatibility with no toxic effect on cells. Easy handling and easy to apply. Ideal for infected root canal and vital pulpotomies in deciduous teeth. It can also be used in conjunctions with gutta percha points and regular root canal sealants.

Shelf life and storage:

- Best before 24 months under proper storage.
- Always close the cap immediately after use.
- Store at a temperature between 10°C to 30°C.