

Calx



Known as: Calcium Hydroxide Paste (Radiolucent)

Description: It is premixed high pH Calcium Hydroxide Water based paste with stabilizer. To be used as a root canal dressing and filling material during endodontic treatment of teeth and as a temporary filling material for pulp capping and pulpotomies.

Packing: 1.5 g x 2 Syr + 3 Applicator Tips

Product Code: 9104/1.5x2

Indication:

- Flowing Canal
- Hard Tissue Formation
- Root Resorption
- Radiolucent

Properties & Advantages

- Ready to use, radiolucent, premixed, non-separating Calcium Hydroxide Paste which is useful in periapical infection.
- Used as a liner.
- It reduces the hypersensitivity of the dentine and helps the formation of secondary dentine.
- Non-Separating Homogenous Mix.
- Ideal Viscosity & Micronized particles of Calcium Hydroxide for Effective Application.
- Easy Insertion in Canal Thanks to “**MAARC ADHERE TECH**”

Shelf life and storage:

- Best before 24 months under proper storage.
- Always close the cap immediately after use.
- Store at a temperature between 10°C to 30°C.