

Endo Paste



Known as: Chelating Material

Description: It is water soluble Ethylene Diamine Tetra Acid (17% EDTA with Carbamide Peroxide) Paste for root canal cleaning.

Packing: 3 g x 1 Syr + 3 Applicator tips

Product Code: 9203/003

Benefits:

- **GLIDE YOUR FILES:** with unmatched lubrication.
- **EASYCOAT:** Coats the file with ease.
- **REMOVE SMEAR LAYER EFFECTIVELY.**
- **EFFERVESCENCE:** Reacts with sodium hypo to flush debris.
- **HYDROPHILIC:** easy wash off.
- **NO WASTAGE:** better stability nondrying & no phase separation.

Indication of Use

- Chelating (Calcium Binding) agent.
- For chemical and mechanical cleaning of the root canal.
- Shaping of root canal.
- To be used with sodium hypochlorite. Do not use on soft tissues light skin.

Storage & Shelf Life:

- 24 months from date of manufacturing.
- Always close the cap immediately after use.
- Store at a temperature between 10°C to 30°C.