

# GSOL



**Known as:** Gutta Percha Solvent

**Description:** D-limonene based G. P solvent for removing Gutta Percha at the time of root canal e-treatment. For use as during re-root canal treatment of teeth to dissolve/soften the gutta-percha cones aiding in the removal from the root canal.

**Packing:** 15 ml

**Product Code:** 3006/015

**Indications:**

- To remove all types of root canal filling materials including Zinc-Oxide-Eugenol cement
- Fitting of the primary gutta percha point in root canal filling techniques using the lateral condensation method.
- Retreatment after root canal filling by gutta percha condensation.
- For removing gutta percha at the time of root canal re-treatment.

**Benefits:**

- G-Sol instantly softens GP (***Rapid Action***)
- Easy removal of Sealer once GP soften
- Orange oil based (***Pleasant flavour***)

**Product Benefits:**

- Non-hazardous/non carcinogenic unlike chloroform/trichloroethylene.
- Less volatile as compared to chloroform / trichloroethylene.
- Based on natural substances d-limonene, does not cause irritation to oral mucosa and periapical tissue.
- Evaporation of solvent is less due to low volatility.

**Shelf life and storage:**

- Best before 24 months under proper storage.
- Always close the cap immediately after use
- Store at a temperature between 10°C to 30°C.



**Manufactured by:**  
**Shiva Products**

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),  
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com