

# T-Fill



**Known as:** Temporary Filling Material

**Description:** T-Fill is ready to use temporary filling material self-curing under humidity for temporary filling of cavities.

**Packing:** 40 g

**Product Code:** 5507/040

**Attributes:**

**Flavor:** Mint

**Colour:** White

**Indication of use:**

- For use as Temporary filling material after Endodontic treatment
- For use as Temporary Inlay material

**Properties:**

- Eugenol Free offer nonirritating properties.
- Higher adhesion to teeth structure.
- Micronized and homogeneously mixed for grain free consistency.
- Gives best adaptability.
- Quick setting in moist conditions
- There is minor but even expansion on setting, so no cracks develop, giving best sealing hence protecting to Pulp.

**Benefits:**

- Excellent Adaptability
- Non-Sticky to Instrument
- No Cracks: Even Expansion Due to Even Grain Size
- Eugenol Free
- Pleasant Flavor: Mint
- Neutral Taste
- Minimum Moisture Contamination: Dual Cap for Air Sealing

**Shelf life and storage:**

- 24 months under proper storage.
- After opening store in refrigerator (do not Freeze).
- Store at a temperature between 10°C to 30°C.



**Manufactured by:**

**Shiva Products**

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),  
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com