

## Beading Wax



**Description:** It is a type of processing wax. Beading Wax is a soft wax ideal for dentists to use to bead round the edge of a functional impression prior to casting.

**Packing:** Pack of 12 strips. It comes in the pack of 12 strips wrapped in non-adhesive Silicone Paper to facilitate easy removal.

**Product Code:** 2309/012

**Thickness:** 6 mm x 6 mm x 250 mm approx.

### Properties & Advantages:

- It is sticky enough to easily adapt at the periphery or edge of alginate impression / final impression by hand pressure.
- It creates the base of the cast and preserves the outer portion of the cast.

### Shelf life and storage:

- Best before 60 months under proper storage.
- Store at Temp between 10°C to 30°C