

# Dental Stone Plaster (Green, Yellow & Off White)



For fabricating custom Trays, mouth guards and vaccum formed splints. For denture repairs. For studycasts. For lab procedures such as articulator cast mounting. For processing full or partial acrylic dentures.

**Packing: (Yellow)** 30kg polyester gunny / paper bag.

**Product Code:** 4001/030

**Packing: (Green):** 30kg polyester gunny / paper bag.

**Product Code:** 4002/030

**Packing: (Off White)** 30kg polyester gunny / paper bag.

**Product Code:** 4005/030

## Type: III - Specification:

- Water/Powder: 30ml/100gm
- Mixing Time: 60 seconds wet
- Working Time: 6 minutes dry
- Setting Time: 8 to 12 minutes
  
- Hard Set (Initial): 30 Minutes
- Wet Compressive Strength: 250 kg/cm<sup>2</sup>
- Dry Compressive Strength: 500 kg/cm<sup>2</sup>
- Setting Expansion (Max): 0.08%

## Shelf life and storage:

Best before 24 months when stored in cool and dry place.