

Green Tracing Sticks



Thickness: 113 x 9 mm dia weighing appx. 10 g

Packing: 10 Sticks

Product Code: 1501/010

Packing: 15 Sticks

Product Code: 1501/015

Packing: 100 Sticks

Product Code: 1501/100

Indication:

- Low fusing material primarily used for border extensions on impression trays.
- High viscosity also makes it ideal when recording the full depth of the sulcus.
- Modify dentate tray by building up material where tooth has been extracted.
- Make functional and relining impressions.

Properties & Advantages:

- Suitable for muscle trimming work.
- Fast and uniform plasticity.
- Quick hardening in cold water.
- Wide working area.

Shelf life and storage:

Best before 60 months under proper storage