

Model Stone Plaster (Blue)



Used making crown and bridge, dies and cast partials, dental arches and opposing models. For making study models, in use where economical yet tough, non-brittle and durable models are required.

Packing: 30kg polyester gunny/paper bag.

Product Code: 4003/030

Type: III - Specification:

- Hard Set (Initial): 30 Minutes
- Setting Expansion: 0.2%
- Compressive 1h. Mpa: 35
- Compressive 24h. Mpa: 75

- Water/ Powder: 26-28ml/100gm
- Mixing Time: 60 seconds wet
- Working Time: 5 minutes dry
- Setting Time: 8-10 minutes

Shelf life and storage:

Best before 24 months when stored in cool and dry place.