

Ortho White Stone Plaster



For crown and bridge, dies and cast partials. For basing dental arches and opposing models.
For making study models.
In use where economical yet tough, non brittle and durable models are required.

Packing: (White) 30kg polyester gunny / paper bag.

Product Code: 4102 / 030

Type: IV – Specification

- Water/ Powder: 24-26ml/100gm
- Mixing Time: 45 seconds wet
- Working Time: 5 minutes dry
- Setting Time: 8-12 minutes

- Hard Set (Initial): 30 Minutes
- Wet Compressive Strength: 400 kg/cm²
- Dry Compressive Strength: 800 kg/cm²
- Linear Expansion (Max): 0.08%

Shelf life and storage:

Best before 24 months when stored in cool and dry place.