

Cold Cure / Rapid Repair Powder & Liquid



Description: Self-cure monomer and polymer can be colorless or pigmented. Self-cure acrylic resin is a chemically setting acrylic resin used in a variety of applications in dentistry in either dough form or brush-on technique or injection moulding technique.

Packing	Product Code	
Cold Cure Powder Clear	6201/100 - 100 g	6201/400 - 400 g
Cold Cure Powder Pink	6202/100 - 100 g	6202/400 - 400 g
Cold Cure Powder Light Veined	6203/100 - 100 g	6203/400 - 400 g
Cold Cure Liquid	6301/100 - 100 ml	6301/400 - 400 ml
Cold Cure Universal Pack Clear	6504/001 - 400 g + 400 ml	
Cold Cure Universal Pack Pink	6505/001 - 400 g + 400 ml	
Cold Cure Universal Pack Light Veined	6506/001 - 400 g + 400 ml	

Indication of use:

- Repair & relining of dentures
- Custom impression trays
- Orthodontic appliances
- Maxillofacial prosthesis
- Provisional inlays and crowns
- Temporary bridges

Shelf life and storage:

- Best before 24 months under proper storage.
- Keep the powder and liquid containers tightly closed to avoid contamination.
- Store at a temperature between 10°C to 30°C.