

# Die Stone (Pink)



For making crown and bridge, dies and cast partials. For basing dental arches and opposing models. For making study models. In use where economical yet tough, non-brittle and durable models are required.

Packing	Product Code
30 kg	4101/030

## Type: IV - Specification:

- Water / Powder 24-26 ml / 100 gm
- Mixing Time 45 seconds wet
- Working Time 5 minutes dry
- Setting Time 8-12 minutes
- Hard Set (Initial) 30 Minutes
- Wet Compressive Strength 400 kg / cm<sup>2</sup>
- Dry Compressive Strength 800 kg / cm<sup>2</sup>
- Linear Expansion (Max) 0.08%

## Shelf life and storage:

- Best before 24 months when stored in cool and dry place.
- Store at Temp between 10°C to 30°C.