

Heat Cure / Cross Linked Powder & Liquid



Description: Heat-cured acrylic denture resin used in conjunction with heat cure acrylic resin powder to make various removable oral & maxillofacial prostheses and pan- facial prosthesis.

| Packing | Product Code | |
|---------------------------------------|-------------------|-------------------|
| Heat Cure Powder Clear | 6001/100 - 100 g | 6001/450 - 450 g |
| Heat Cure Powder Pink | 6002/100 - 100 g | 6002/450 - 450 g |
| Heat Cure Powder Light Veined | 6003/100 - 100 g | 6003/450 - 450 g |
| Heat Cure Liquid | 6101/100 - 100 ml | 6101/225 - 225 ml |
| Heat Cure Universal Pack Clear | 6501/001 | |
| Heat Cure Universal Pack Pink | 6502/001 | |
| Heat Cure Universal Pack Light Veined | 6503/001 | |

Properties & Advantages:

- It is a 2-part rigid dental acrylic used for making dental appliances.
- Contains both powder and liquid components
- Used for full or partial, upper or lower dentures
- Simple mixing, easy to process.
- Cross Linking ensures superior quality.

Shelf life and storage:

- Best before 36 months under proper storage.
- Store at a temperature between 10°C to 30°C.