

# Mock-Up Wax



Packing	Product Code
12 Sheets	2201/225
200 g	2201/200

**Thickness:** 1.3 mm - 1.5 mm (Custom Size / Thickness Catered)

**Size:** 160 mm x 90 mm

**Packing:** Each Packet of 'MAARC' Mockup Wax contains 12 sheets wrapped in non-adhesive Silicone Paper to facilitate easy removal.

**Color:** White

## Product features:

- It simulates the final outcome. The result is temporary and can be changed according to the patient's desire.
- It can be repeated as many times as needed in order to select the best possible solution
- It helps create trust between the patient and the dentist.
- It allows the dentist to evaluate the patient's desires and change the treatment plan accordingly.

## Shelf life and storage:

- Store at Temp between 10°C to 40°C.
- Best before 60 months under proper storage.