

Soft-Liner Kit



Description: Soft denture reline material

Packing	Product Code
50 g + 50 ml	8005/001

Intended use / Indications:

- For the temporary re-basing (relining) of acrylic dentures.
- Conditioning Tissue.

Properties & Advantages:

- Comfortable Tissue Conditioning
- Absorbs Shock Efficiently
- Enhances Denture Stability
- Excellent Denture Adhesion

Shelf life and storage:

- Best before 24 months under proper storage.
- Store at a temperature between 4°C to 25°C.