

Arch Bar



Description: It is made of 316 grade SS suitable for rigid fixation, techniques in the dentate patient begin with fixation of the occlusion. This ensures that the patients maintain their preoperative occlusal status. There are several techniques to providing mandibulomaxillary fixation (MMF). Many surgeons agree that the gold standards in MMF is, the use of arch bars

Packing	Packing code
Arch bar 12 inch long	5508/012

Intended use / Indications:

- For temporary fragment stabilization in emergency cases before definitive treatment. As a tension band in combination with rigid internal fixation
- For long-term fixation in conservative treatment
- For fixation of avulsed teeth and alveolar crest fractures