

Calcium Hydroxide Powder



Known as: Pulp Capping Material

Description: Calcium hydroxide is used as an agent to induce apexification in pediatric endodontics. The product may also be used as a pulp capping material and as a base under non-adhesive restorations.

Packing	Product Code
15g	5002/015
100gm	5002/100

Intended use / Indications:

- Ideal for temporary root filling
- Pulp capping
- Pulp protection
- Provides insulation in deep cavities

Properties & Advantages:

- Provides bactericidal effect on carious dentine.
- Paste delivered to root apex provides pulp insulation.
- Protects the vital pulp from infections and further decay.
- Can be easily removed from root canal by simple irrigation.

Shelf life and storage:

- Best before 24 months under proper storage.
- Store at room temperature between 10° C to 30° C.