

Carding Wax / Teeth Wax



Packing: 100 g

Product code: 2804/100

We manufacture two types of carding wax for Pouring and for Extrusion.

Properties & Advantages:

- Both the waxes are very stable at temperature from 12°C to 40°C.
- Pouring waxes are suitable for automatic machines for tray filling.
- Wax firmly holds the teeth and while teeth is pulled negligible amount of wax comes along with it.
- Teeth can be re-adhere gaining similar strength

Shelf life and storage:

- Best before 60 months under proper storage.