

CeraSeal-B



Known as: Bio Ceramic Sealer

Description: MAARC CeraSeal-B is MTA based Bio-ceramic Root Canal Sealer.

Packing	Product Code
1 g powder, 5 ml liquid, measuring scoop	8501/001

Intended Use/ Indications: Permanent obturation of root canal of teeth with obturation point.

Properties & Advantages:

- It is resin and monomer free and ensures Zero shrinkage and biocompatibility.
- When mixed, produces mixture with desired flow and adhesion.
- Its high pH and Bioactive properties help remineralization due to hydroxy-apatite formation.
- It has good radiopacity & Due to its purity it will not stain teeth.

Shelf life and storage:

- Best before 24 months under proper storage
- Store at a temperature between 10°C to 30°C.