


# Die Spacer



**Description:** Die Spacer is used to provide sufficient space between the die and crown to accommodate cements used during seating.

| Packing   | Type   | Product Code |
|---|--------|--------------|
| 100ml glass bottle with 10ml applicator bottle+50ml Thinner | Golden | 8101/100     |
| 100ml glass bottle with 10ml applicator bottle+50ml Thinner | Silver | 8102/100     |

**Colour:** Silver/Golden

### Intended Use/ Indication:

Die Spacer is a coloured, paint-on material applied to dies that adds a desired dimension to the surface of the die to function as a “spacer.”

### Properties and Advantages:

- It produces a homogenous coating of 5 ml in single application on the die for crown and bridge restoration thereby assuring accurate seating of restorations.
- The film does not disassociate from the die while building wax pattern due to heating.
- It is rapid setting and thereby reduces working time.

### Shelf life and Storage:

- 24 months under proper storage
- Always close the cap immediately after use.



Manufactured by:

**Shiva Products**

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),  
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com