

# Eugenol



Known as: Eugenol based Temporary Filling material.

Description: Eugenol, also called clove oil, is an aromatic oil extracted from cloves

Packing	Product Code
15 ml	3002/015
30 ml	3002/030
110 ml	3002/110

## Intended use / Indications

- For temporary relief of localized dental pain when combine with zinc oxide to form a sedative dressing.
- Chemically pure with increased effectiveness.
- Forms a uniform paste consistency when mixed with zinc oxide.
- Forms a sedative dressing or lining under non composite resin when mixed with zinc oxide

## Shelf life and storage:

- Best before 24 months under proper storage
- Store at room temperature between 10° C to 30° C.