

# Find Fast



**Known as:** Caries Indicator

## Description:

Find Fast is an important tool for the detection of decay /caries through the method of pigmentation of dental caries. Find Fast clearly differentiates clinically the two layers of carious & normal dentin by staining the outer carious dentin scarlet red yet not staining the inner normal dentin. The distinct staining of only the outer, infected carious layer provides a visible, clinically precise guide for caries removal.

| Packing | Product Code |
|---------|--------------|
| 15 ml   | 3005/015     |

## Intended use / Indications:

Detecting carious dentine during cavity preparation in restorative dentistry. Detection of dental caries in places difficult to identify and prevent excessive removal of normal dentin.

## Properties & Advantages:

Identifies carious dentin ensures the maximum amount of inner dentin is preserved and infected dentin is removed.

Easily locate root canal orifices and calcified canals.

## Shelf life and storage:

- Best before 24 months under proper storage.