

GSOL



Known as: Gutta Percha Solvent

Description: It is an orange oil based Gutta Percha solvent used for removing gutta percha at the time of Re-root canal treatment. It is non-hazardous as it is extracted from naturally existing substances which are considered as non-carcinogenic. Orange oil does not irritate soft tissue or periapical tissues and also offers pleasant fragrance. It does not evaporate as rapidly as chloroform or other solvents.

Packing	Product Code
15 ml	3006/015

Intended use/ Indications:

For use as during re-root canal treatment of teeth to dissolve/soften the gutta-percha cones aiding in the removal from the root canal. Re-root Canal treatment of teeth.

Properties & Advantages:

- Extracted from natural resources.
- Softens the GP Points aiding in their retrieval by mechanical methods.

Shelf life and storage:

- Best before 24 months under proper storage.
- Always close the cap immediately after use
- Store at a temperature between 10°C to 30°C.