

## Model Stone Plaster (Blue)



**Description:** For crown and bridge, dies and cast partials. For basing dental arches and opposing models. For making study models. In use where economical yet tough, non-brittle and durable models are required

Packing	Product Code
(White) 30kg polyester gunny / paper bag.	4003/030

### Type: III – Specification:

- Hard Set (Initial): 30 Minutes
- Setting Expansion: 0.2%
- Compressive 1h. Mpa: 35
- Compressive 24h. Mpa: 75
- Water/ Powder: 26-28ml/100gm
- Mixing Time: 60 seconds wet
- Working Time: 5 minutes dry
- Setting Time: 8-10 minutes
- 

### Shelf Life and Storage:

- Best before 24 months under proper storage.
- Store at a temperature between 10°C to 30°C.