

Ni – Cr Alloy / White Metal for Cast Restorations & Crown & Bridge (Non-Precious Alloy)



Packing	Product Code
500 gms	7001/500

Composition: Ni - 68%, Cr - 15%, Cu - 15% and Rest - 2%

Technical Data: Coefficient of expansion (0-500°C): range from 13.9-14.1µm/m°C

Properties & Advantages:

- It is beryllium free alloy with excellent castability.
- It is reusable with addition of appx 50% fresh metal during recasting.
- It can be cast in open flame.
- It has got good burnish ability and is easy to finish.
- It can be used for making cast crown bridge and pontics and frames for hard resin crowns.



Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com