

Sealmax-R



Known as: Resin based root canal sealer

Description: Sealmax-R is a resin-based root canal sealing material supplied in the form of powder and liquid. The powder and liquid are mixed to produce a root canal sealing material of excellent sealing capacity and tissue compatibility. As a result of its close adaptation to the walls and hydrophobic nature. Sealmax-R allows tight sealing of root canal during obturation.

Packing	Product Code
10 g powder, 10 ml liquid, measuring scoop and mixing pad	8001/001

Intended Use/ Indications: Sealing of root canals during obturation in the secondary dentition with or without the aid of obturation points (gutta-percha).

Properties & Advantages

- Eugenol Free
- Insoluble in water and non-degradable.
- Minimum shrinkage at room temperature.
- Grain Free
- Hydrophobic
- Bacteriostatic
- Negligible Shrinkage
- Complete Impervious Seal
- Excellent Adhesion
- Tissue Tolerant, Most Ideal Viscosity
- Best Intracanal Sealing
- Excellent Radio Opacity

Shelf life and storage:

- 24 months under proper storage.
- Always close the cap immediately after use.
- Store at a temperature between 10°C to 30°C.