

Twin Kleen [The All-in-One Irrigating Solution]



Known as: Single step – Chelation and Irrigation

Description: Twin Kleen is a Mild Chelating agent HEBP with short term compatibility with sodium hypochlorite. The combination eliminates the debris impaction in the anatomical irregularity. The combination also reduces the AHTD (Accumulated Hard Tissue Debris) and prevents smear layer formation during rotary instrumentation. Mixing Tween Kleen with sodium hypochlorite. One can avoid the irrigation sequence and the combination can be used as one irrigating solution. The combination can be used during instrumentation and as a final rinse solution.

Packing	Product Code
100 Capsules	3014/100
25 Capsules	3014/025
50 Capsules	3014/050

Indications:

Twin Kleen is an additive to endodontic irrigating solutions for all endodontic procedure during root canal instrumentation and as a final rinse solution.

Properties & Advantages

- Antibacterial solution
- Mild Chelating Agent
- Proteolytic Agent
- Lubricating Agent
- Prevents smear layer

Shelf life and storage:

- Best before 24 months under proper storage.
- Store at a temperature between 10°C to 30°C.



Manufactured by:
Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com