

Zirprophy



Known as: Prophylaxis Polishing Paste

Contains: Zirconium Silicate

Packing	Product Code
45 g	9001/045
340 g	9001/340

Properties & Advantages:

- Zirprophy contains Zirconium Silicate.
- It does not break down as pumice from coarse to fine but rather removes stain smoothly due to fine particles.
- Zirprophy is to be used light handedly to avoid damage to enamel
- Thixotropic Oil Free
- Splatter Free
- Fluoride Free
- Spearmint Flavor
- Zirconium silicate abrasive media
- Easy removal of smoking and tartar stains
- Controlled particle size (Minimum post scaling sensitivity)
- Better polish
- Fluoride Free

Shelf life and storage:

- Best before 24 months under proper storage.
- Always close the cap immediately after use.
- Do not freeze.