

Beading Wax



Description: It is a type of processing wax. Beading Wax is a soft wax ideal for dentists to use to bead round the edge of a functional impression prior to casting.

Packing: It comes in the pack of 12 strips wrapped in non-adhesive Silicone Paper to facilitate easy removal.

Product Code	Packing
2309/012	Pack of 12 strips

Thickness: 6 mm x 6 mm x 250 mm approx.

Intended use / Indications: Beading Wax is adapted around the impression borders to create the land area of the cast.

Properties & Advantages:

- It is sticky enough to easily adapt at the periphery or edge of alginate impression / final impression by hand pressure.
- It creates the base of the cast and preserves the outer portion of the cast.

Shelf life and storage:

- Store at temp. between 10°C to 30°C
- Best before 60 months under proper storage.