

Carving Wax / Wax Block



'MAARC' Carving wax comes in small blocks slugs and in various colors like Translucent White, Yellow, Red, Blue.

Description

Carving Waxes are made from high-quality ingredients and widely used in dentistry or dental lab.

Packing	Product Code
40 Blocks - White	2501/040
40 Blocks - Red	2502/040
40 Blocks - Blue	2503/040
40 Blocks - Yellow	2504/040

Thickness: 38 x 13 x 13 mm

Intended use / Indications:

It is applied to crown and bridges and is appreciated for its tough, hard, and bubble-free attributes as well as it does not chip or flake.

Properties & Advantages:

- The Wax Blocks are easy to carve.
- The special feature being the wax is sufficiently hard that it does not chip off at edges on applying pressure.
- Free from air bubbles.

Shelf life and storage:

- Store at Temp between 10°C to 30°C
- Best before 60 months under proper storage.



Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com