

Clasp wax



Also known as: Cast Partial Waxes

Description:

Clasp Profile waxes are primarily used in the wax-up technique; this can include sculpturing, build up and inlay work.

Packing	Product Code			
10 Sheets	2901/010	2904/010	2907/010	2910/010
	2902/010	2905/010	2908/010	
	2903/010	2906/010	2909/010	

Intended use / Indications:

Wax clasp profiles for molars and premolars are designed to ensure good hold and stability for dentures. Half-teardrop shaped cross-sections are effective at preventing food residues from getting stuck. All clasp profiles have an optimal anatomical shape and are made of easily adapted wax. Adhesion to the investment material model is optimal.

Packing: 10 Pcs

Shelf life and storage:

- Store at Temp between 10°C to 30°C
- Best before 60 months under proper storage.



Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com