

Inlay Sticks



Description:

Inlay/Onlay wax are type-1 waxes which soften at a temperature that is not hazardous to the dental pulp tissue, hence can be used to prepare patterns of indirect-restorations in the patient's mouth by direct-wax technique. They have uniform workability over a wide temperature range and facilitate accurate adaptation to the prepared cavity.

Thickness: L 114 x 8 mm dia of Big Sticks

Packing	Product Code
10 Sticks (70 g)	2403/010

Intended Use/ Indications:

Used for the construction of wax pattern using direct wax technique, for fabrication of inlay & onlay restorations, fabrication of onlays and full coverage crowns by indirect wax technique on the die models.

Shelf life and storage:

- Store at Temp between 10°C to 30°C
- Best before 60 months under proper storage



Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarc dental.com