

Spacer Wax



Thickness: 0.4 mm (+/-0.1 mm)

Packing: Upper and Lower Spacer Wax are packed separately. Wrapped in non-adhesive Silicone Paper to facilitate easy removal.

Packing	Product Code
Upper-25 Sheets	2601/025
Lower-25 Sheets	2602/025

Intended use / Indications: To provide space for final impression material

Properties & Advantages:

- Thin wax sheets of 0.3mm thickness - shape of Upper and Lower Base Plates
- It increases the thickness of wax rims while making dentures using base plates.
- Provide uniform thickness of impression material.
- Use of spacer wax prevents direct contact of acrylic trays with soft tissues

Shelf life and storage:

- Store at Temp between 10°C to 30°C
- Best before 60 months under proper storage.



Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com