

Utility Wax



Description:

Utility Wax are in form of ropes / wires, easily shaped waxes that are especially suited for fill-ins and wax pattern repair. It can also be used to easily correct mistakes by “blending-in” with your original wax part.

Thickness: 5.0 diameter round x 30cm long

Packing: 500gm Box / 55 Rods

Product Code	Packing
2308/012	12 Strips
2308/250	250 g

Intended use / Indications:

Utility wax is used for protecting and preserving the periphery edge in removable prosthetics. The application of utility wax on the impression tray allows a perfect transfer of the periphery edge from the impression onto the model.

Properties & Advantages:

1. Extremely pliable at room temperature.
2. It can also be used to fill sprue bases.
3. The wax can be adapted at room temperature and can be used without pre-heating.

Shelf life and storage:

- Store at Temp between 10°C to 30°C
- Best before 60 months under proper storage.