

# GSOL



**Known as:** Gutta Percha Solvent

Packing	Product Code
15 ml	3006/015

### Intended use/ Indications:

GSOL is intended for use during repeated root canal treatments to remove gutta-percha cones from the root canal.

### Product Description & Properties:

It is an orange oil (D-limonene) based Gutta Percha Solvent used for removing gutta percha at the time of Re-root canal treatment.

### Shelf life and storage:

- Best before 24 months from the date of manufacturing under proper storage.
- Always close the cap immediately after use to prevent contamination.
- Store at a temperature between 10°C to 30°C.