

Articulating Paper



Description: Articulating paper is commonly used in dentistry to identify contactpoints between the maxillary and mandibular teeth during all forms of natural tooth occlusal adjustments and dental prosthesis insertions.

Packing	Product Code
40 μ- Red - 1 Box containing 10 books	5515/001
40 μ- Blue - 1 Box containing 10 books	5519/001
40 μ- Red / Blue - 1 Box containing 10 books	5511/001
100 μ- Red - 1 Box containing 10 books	5517/001
100 μ- Blue - 1 Box containing 10 books	5521/001
100 μ- Red / Blue - 1 Box containing 10 books	5513/001

Intended use / Indications: Used as diagnostic device in dentistry to highlight occlusal contacts and the distribution of occlusal forces.

Properties & Advantages:

- Maarc articulating papers are coated with easy transferable ink, so it accurately highlights occlusal contacts.
- It comes in various thickness.
- Maarc articulating paper comes with progressive color transfer and aids in visualizing static occlusion.
- Because of its low thickness and flexibility, it adapts perfectly to cusps and fossae, making it extremely precise and therefore suitable for more delicate preparations.

Thickness: 40μ / 100μ

Shelf life and storage:

- Store at Temp between 10°C to 40°C.
- Best before 60 months under proper storage.



Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com