

Boxing Wax



Description:

Used for boxing impressions prior to pouring the model. A tacky wax that results in perfect adhesion without the application of heat.

Thickness: 160 x 45 x 1.5 mm

Packing: Pack of 6 Sheets wrapped in silicon paper to facilitate easy removal.

Packing	Product Code
6 inch-6 Sheets	2203/006
12 inch-6 Sheets	2204/006

Intended use/Indications:

Boxing Wax is used to build up vertical walls around the impression in order to pour the gypsum product to make a cast base.

Advantages:

- It is soft and tacky for boxing in impressions.
- Adheres easily to most of the surfaces.

Shelf life and storage:

- Store at Temp between 10°C to 40°C
- Best before 60 months under proper storage



Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com