

Casting Wax Sheets



Description: Dental casting waxes are primarily used in the wax-up technique; this can include sculpturing, build up and inlay work.

Thickness: 0.3 to 0.6 mm

Size: 160 mm x 90 mm

Colour: 'Maarc' Casting Wax sheets are thin wax sheets. Available in Green color for cobalt chrome cast partial dentures.

Packing: Available in 200 g pack. The Sheets are wrapped in non-adhesive Silicone Paper to facilitate easy removal.

Packing	Product Code
200 g	2206/200

Intended use / Indications: Casting wax is used to form the pattern for a partial denture framework in a dental casting alloy such as Cobalt-Chromium.

Properties & Advantages:

- The Sheets are sufficiently thin to have thinner casting as possible without defects to reduce machining.
- They easily burn out during castings.

Shelf life and storage:

- Store at Temp between 10°C to 30°C.
- Best before 60 months under proper storage