

Cold Cure / Rapid Repair Powder & Liquid



Product Type	Packing	Product Code
Cold Cure Powder Clear	100 gm, 400 gm, 500 gm, 1 kg, 3 kg, 5 kg	6201/100, 6201/400, 6201/500, 6201/001, 6201/003, 6201/005
Cold Cure Powder Pink	100 gm, 400 gm, 500 gm, 1 kg, 3 kg, 5 kg	6202/100, 6202/400, 6202/500, 6202/001, 6202/003, 6202/005
Cold Cure Powder Light Veined	100 gm, 400 gm, 500 gm, 1 kg, 3 kg, 5 kg	6203/100, 6203/400, 6203/500, 6203/001, 6203/003, 6203/005
Cold Cure Powder Dark Veined	100 gm, 400 gm, 500 gm, 1 kg, 3 kg, 5 kg	6204/100, 6204/400, 6204/500, 6204/001, 6204/003, 6204/005
Cold Cure Liquid	100 ml, 400 ml, 500 ml, 1 ltr, 3 ltr	6301/100, 6301/400, 6301/500, 6301/001, 6301/003
Cold Cure Universal Pack Clear	400 g + 400 ml	6504/001
Cold Cure Universal Pack Pink	400 g + 400 ml	6505/001
Cold Cure Universal Pack Light Veined	400 g + 400 ml	6506/001

Intended Use /Indication:

- Repair & relining of dentures
- Custom impression trays
- Orthodontic appliances
- Maxillofacial prosthesis
- Provisional inlays and crowns
- Temporary bridges.

Shelf Life and Storage:

- Best before 36 months under proper storage.
- Store at a temperature between 10°C to 30°C.
- Keep the powder and liquid containers tightly closed to avoid contamination.



Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com