

# Cold Mould Seal / Separating Medium



**Description:** Cold Mould Seal is a specially developed formulation of separating solution for universal use, suitable for both packing / pressing technique & cast moulding.

Packing	Product Code
110 ml	6401/110
5 liters	6401/005

## Intended use / Indications:

Used when great accuracy is not required eg.

- Flasking of denture
- Articulating the model
- Separation of model from base during repair works

## Properties & Advantages:

- Makes tough film with elasticity which is unbreakable in pressure.
- Film does not chip off from plaster at pressure. High quality of the raw materials guarantees optimum separation with both techniques

## Shelf life and storage:

- Best before 36 months under proper storage.
- Store at a temperature between 10°C to 30°C.