

## ETC-37



**Known as:** Etching Materials (37% Phosphoric Acid)

**Description:** ETC - 37 is a specially formulated surface etchant which can be used for preparing the tooth enamel and dentin to receive bonded restorations. Supplied in the form of a thixotropic gel/liquid it provides excellent control over the flow of the etchant. The device allows easy manipulation by the dentist and etch the surface to precise extent. It can also be used to etch the surfaces of enamel prior to restorations with resin based restorative materials.

Packing	Product codes
3 g x 1 Syr + 3 Applicator tips	9401/003
25 g x 1 Syr + 10 Applicator tips + 1 Connector	9401/025
5 g x 3 Syr + 10 Applicator tips	9401/5x3
15 ml x 1 liquid bottle refill	9402/015

### Intended use / Indications:

1. Surface etching of the surfaces of teeth which are to receive bonded restorations.
2. Surface etching of the prepared access cavity in endodontically treated teeth to receive bonded restorations.
3. Etching of both enamel and dentinal surfaces of teeth prior to resin based restorations.
4. Etching of the surfaces of teeth prior to bonding in orthodontic treatment of teeth.

### Shelf life and storage:

- Best before 24 months under proper storage.
- Store at a temperature between 10°C to 30°C.



Manufactured by:  
**Shiva Products**

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),  
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com