

# Endo Paste



**Known as:** Chelating Material

**Description:** It is water soluble Ethylene Diamine Tetra Acid (17% EDTA with Carbamide Peroxide) Paste for root canal cleaning.

Packing	Product Code
3 g x 1 Syr + 3 Applicator tips	9203/003

### Properties & Advantages:

1. Effective chelating, lubrication, deburring water-soluble material
2. It facilitates root canal preparation

### Intended Use/Indication:

For use as root canal lubricant during mechanical instrumentation in endodontic treatment of teeth. For use as a root canal conditioner during chemico-mechanical preparation of calcified root canal and shaping the canal, aiding in the softening of the calcifications and removal of debris from the prepared tooth.

### Storage & Shelf Life:

- 24 months from date of manufacturing.
- Always close the cap immediately after use.
- Store at a temperature between 10°C to 30°C.