

Endogel-RC & Endo-L



Known as: Chelating Material

Description:

Endogel-RC: Endogel-RC is a ready to use water-soluble Ethylene Diamine Tetracetic Acid (17% EDTA) in Gel form for root canal cleaning, packed in multi-dose plastic dispensing syringe. It is both lubricating and chelating material in gel form, to be used during cleaning and shaping of the root canals during endodontic treatment. This is non sterile device & does not require sterilization while being used.

Endo-L: Endo-L is a ready to use water-soluble Ethylene Diamine Tetracetic Acid (17% EDTA) in Liquid form for root canal cleaning, packed in multi-dose plastic dispensing bottle. It is both lubricating and chelating material in liquid form, to be used during cleaning and shaping of the root canals during endodontic treatment. This is non sterile device & does not require sterilization while being used.

Packing	Product Code
Endogel-RC: 3 g x 1 syr + 3 Applicator tips	9202/003
Endo-L: 100ml Bottle Refill	9201/100

Intended use/ Indications:

For use as root canal lubricant during mechanical instrumentation in endodontic treatment of teeth. For use as a root canal conditioner during chemico-mechanical preparation of calcified root canal and shaping the canal, aiding in the softening of the calcifications and removal of debris from the prepared tooth.

Properties and Advantages:

- Effective chelating.
- Lubrication.
- Deburring water-soluble material.
- 17% Liquid based EDTA.
- Used to flush out the debris from canal.
- Removes the smear layer from dentinal walls.
- It facilitates root canal preparation.
- Also helps in removal of calcium hydroxide oil based from the canals.
- Mechanically activated using an endo brush.

Shelf Life and Storage:

- Best before 24 months under proper storage.
- Always close the cap immediately after use.



Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com

- Store at a temperature between 10°C to 30°C.



Manufactured by:

Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com