

## Heat Cure / Cross Linked Powder & Liquid



Product Type	Packing	Product Code
Heat Cure Polymer - Clear	100 gm, 400 gm, 500 gm, 1 kg, 3 kg, 5 kg	6001/100, 6001/400, 6001/500, 6001/001, 6001/003, 6001/005
Heat Cure Powder - Pink	100 gm, 400 gm, 500 gm, 1 kg, 3 kg, 5 kg	6002/100, 6002/400, 6002/500, 6002/001, 6002/003, 6002/005
Heat Cure Powder - Light Veined	100 gm, 400 gm, 500 gm, 1 kg, 3 kg, 5 kg	6003/100, 6003/400, 6003/500, 6003/001, 6003/003, 6003/005
Heat Cure Powder - Dark Veined	100 gm, 400 gm, 500 gm, 1 kg, 3 kg, 5 kg	6004/100, 6004/400, 6004/500, 6004/001, 6004/003, 6004/005
Heat Cure Liquid	100 ml, 400 ml, 500 ml, 1 ltr, 3 ltr	6101/100, 6101/400, 6101/500, 6101/001, 6101/003
Heat Cure Universal Pack Clear	450 gms + 225 ml + 110 ml CMSL	6501/001
Heat Cure Universal Pack Pink	450 gms + 225 ml + 110 ml CMSL	6502/001
Heat Cure Universal Pack Light Veined	450 gms + 225 ml + 110 ml CMSL	6503/001

### Intended Use/ Indication:

- Removable complete denture
- Removable partial denture
- Orthodontic appliances
- Oral & Maxillofacial prosthesis
- Obturators
- Relining & Rebasing of dentures.

### Shelf Life and Storage:

- Best before 36 months under proper storage.
- Store at a temperature between 10°C to 30°C.
- Keep the powder and liquid containers tightly closed to avoid contamination.



### Manufactured by:

**Shiva Products**

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),  
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com