

Shellac Base Plates



Description: Shellac Material are widely used in prosthetic dentistry, as temporary baseplate for occlusal rims and trial dentures.

Thickness: Upper / Lower Base Plate: 1.3mm, 1.7mm, 2.0mm, 2.5mm and 3.0mm (+/- 0.15mm) (Custom Thickness Catered)

Color: Brown /Yellow/ Pink / Aluminum

Packing: Each Packet of 'MAARC' Base Plate contains 8 upper and 4 Lower Plates, 50/100 Pcs. (Custom Packing Available)

Packing	Product Code
Brown Base Plate (8 Upper + 4 Lower)	1001/084
Yellow base Plate (8 Upper + 4 Lower)	1024/084

Type: It comes in different formulations for different climatic conditions (LMP / HMP / AI. Filled).

Intended use / Indications:

- Individual impression trays
- Making dentures: Base material for bite impressions and tooth lineups

Properties & Advantage:

- They are ready to use Maxillary & Mandibular plates - 8 upper & 4 lower.
- Optimum heat resistance to minimize overheating of the Plate.
- Easily adaptable on heating and easy to work with.
- Retain the shape stability and do not warp back on cooling.
- Optimum mechanical strength to resist breaking during Try-in denture.
- Does not adhere to plaster.

Shelf life and storage:

- Best before 24 months.
- Store below 30°C under proper storage and avoid from direct sunlight.



Manufactured by:
Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com