

# Sticky Wax Sticks



**Description:** Wax to fix removable prosthesis. Its high adherent properties make of this wax an excellent material to fix resin teeth to the anterior walls and to fix these walls to the model.

Packing	Product Code
10 pcs	2701/010
200 g	2701/200

## Intended use / Indications:

- It is used to hold or weld materials or pieces together.
- The Sticky Wax can also be used to connect or repairs wax patterns.

## Properties & Advantages:

- It is blending of various waxes and resins and it is yellow in color.
- It melts when heated and adheres to the surface over which it is applied.
- It holds the components prior to fixation and gives a clean burnout.

## Shelf life and storage:

- Best before 60 months under proper storage.
- Store at a temperature between 10°C to 30°C.