

T-Fill



Known as: Temporary Non-Eugenol Cement (Temporary Filling Material)

Description: T-Fill is ready to use temporary filling material self-curing under humidity for temporary filling of cavities.

Packing	Product Code
40 g	5507/040

Product Attributes:

Flavor: Mint

Colour: White

Intended Use/ Indication:

- For use as Temporary filling material after Endodontic treatment
- For use as Temporary Inlay material

Properties & Advantages:

- Eugenol Free offer nonirritating properties.
- Higher adhesion to teeth structure.
- Micronized and homogeneously mixed for grain free consistency.
- Gives best adaptability.
- Quick setting in moist conditions
- Excellent Adaptability
- Non-Sticky to Instrument
- No Cracks: Even Expansion Due to Even Grain Size
- Pleasant Flavor: Mint
- Neutral Taste
- Minimum Moisture Contamination: Dual Cap for Air Sealing

Shelf life and storage:

- 24 months under proper storage.
- After opening store in refrigerator (do not Freeze).
- Store at a temperature between 10°C to 30°C.