

Inlay/ Onlay Pattern Wax



Description: Inlay/Onlay wax are type-1 waxes which soften at a temperature that is not hazardous to the dental pulp tissue, hence can be used to prepare patterns of indirect-restorations in the patient's mouth by direct-wax technique. They have uniform workability over a wide temperature range and facilitate accurate adaptation to the prepared cavity.

Packaging	Product Code	Description
Cervical Wax (1) Orange (200g)	2401/200	Soft type with degree of fluidity that will allow it to flow and fit perfectly in cervical areas, high adaptation and stability to surfaces. It will wet model surface without shrinkage.
Cervical Wax (2) Red/Blue/Light Blue (200g)	2402/200, 2409/200, 2410/200	Harder than orange wax.
Cervical Wax (3) Blue/Red/Ivory (200g)	2403/200, 2412/200, 2413/200	Standard inlay wax wherein the pattern comes out easily from the plaster surface.
Cervical Wax (4) Green (200g)	2404/200	Harder than red wax used where hard pattern waxes are preferred.
Cervical Wax (5) Violet/Ivory (200g)	2405/200, 2411/200	Harder than green wax – used for copying of patterns to add mechanical strength and as undercut wax.
Cervical Wax (6) Coffee (200g)	2406/200	Used for undercuts and adding mechanical stability to wax patterns of jacket crowns.
Cervical Wax (7) Magenta (200g)	2407/200	Hard inlay wax. With medium flow for built up.
Presentation Orange Inlay Wax (200g)	2408/200	A white, enamel-like wax used to make a more aesthetic case presentation (wax-up). Easy carving properties. Also aids in the fabrication of temporaries.

Properties & Advantages:

- Softer waxes with high flow are more flexible and are often used.

Shelf life and storage:

- Store at Temp between 10°C to 30°C.
- Best before 60 months under proper storage.



Manufactured by:
Shiva Products

HO: 214, Bluebell Indl. Estate, Chinchpada, Vasai (East),
Dist. Palghar, 401 208, Maharashtra, India

+919967806543 | care@shivaproducts.net

www.maarcidental.com